

NAVIGATOR
Dilution Control System

NAVIGATOR
Dilution Control System

PERISEPT

PERISEPT

PERISEPT

62

**SPORICIDAL
DISINFECTANT
CLEANER**

62



62

**SPORICIDAL
DISINFECTANT
CLEANER**

HOSPITAL DISINFECTANT

**SPORICIDAL
DISINFECTANT
CLEANER**

HOSPITAL DISINFECTANT

Before Using This Product,
Please Read This Entire Label Carefully

Before Using This Product,
Please Read This Entire Label Carefully

Active Ingredients:
Hydrogen Peroxide.....27.3%
Peroxyacetic Acid.....5.9%
Inert Ingredients:.....66.8%
Total.....100.0%

KEEP OUT OF REACH OF CHILDREN
See back panel of label for additional
precautionary statements and first aid.
DANGER
PELIGRO
E.P.A. Reg. No. 10324-214-12120
E.P.A. Est. No. 8325-PA-01

Batch Code: 914492 031220 0800

NET CONTENTS: 67.8 fl. oz. (2 Liters)

NAVIGATOR
Dilution Control System

NAVIGATOR
Dilution Control System

PERISEPT

PERISEPT

PERISEPT

62

**SPORICIDAL
DISINFECTANT
CLEANER**

62



62

**SPORICIDAL
DISINFECTANT
CLEANER**

HOSPITAL DISINFECTANT

**SPORICIDAL
DISINFECTANT
CLEANER**

HOSPITAL DISINFECTANT

Before Using This Product,
Please Read This Entire Label Carefully

Before Using This Product,
Please Read This Entire Label Carefully

Active Ingredients:
Hydrogen Peroxide.....27.3%
Peroxyacetic Acid.....5.9%
Inert Ingredients:.....66.8%
Total.....100.0%

KEEP OUT OF REACH OF CHILDREN
See back panel of label for additional
precautionary statements and first aid.
DANGER
PELIGRO
E.P.A. Reg. No. 10324-214-12120
E.P.A. Est. No. 8325-PA-01

Batch Code: 914493 031220 0800

NET CONTENTS: 67.8 fl. oz. (2 Liters)

- Effective in 2 minutes against *C. difficile* spores.
- Cuts cleaning time.
- Is a heavy duty disinfectant cleaner that cleans and disinfects in one labor saving step.
- Leaves no visible residue.
- No rinsing.
- Effective in 10 minutes against *Mycobacterium bovis* (Tb).
- Good for use with microfiber cloths.
- May be used to clean and disinfect finished floors.
- Cleans and disinfects without dulling gloss.
- Economical concentrate.
- Intended for use with the NAVIGATOR Dilution Control System.
- Closed loop automated dispensing reduces employee exposure to concentrate product.
- NAVIGATOR Dilution Control System makes accurate dispensing quick and easy.
- Designed for daily use on common materials found in hospitals.

Daily use product with *Clostridium difficile* spore efficacy allows for product standardization and eliminates need for bleach. Effective against Multidrug Resistant Organisms (MDROs) *Staphylococcus aureus*, Resistant to Methicillin (MRSA), *Staphylococcus aureus*, Community Associated Methicillin Resistant (CA-MRSA), *Staphylococcus aureus*, Intermediate Vancomycin Resistance (VISA), *Enterococcus faecalis* Resistant to Vancomycin (VRE), *Escherichia coli* Extended-Spectrum Beta-Lactamase resistant (ESBL), *Klebsiella pneumoniae* Carbapenem resistant (KPC).

MATERIAL COMPATIBILITY

Not recommended for use on copper, brass, granite, marble or zinc. Do not use on unsealed/uncoated marble or unsealed/uncoated terrazzo floors.

NOTE: This product is compatible with the listed materials. If product is intended to be used on any other surface, it is recommended that you apply product to a smaller test area to determine compatibility before proceeding with its use.

SURFACES: Finished floors, aluminum, chrome, glazed ceramic, laminated surfaces, plastic, painted surfaces, stainless steel, vinyl and plastic upholstery, shower stalls, bathroom fixtures, basins, hospital beds, bed railings, wheelchairs, x-ray tables, medical equipment surfaces and other hard, non-porous surfaces.

For use in Hospitals, Nursing Homes, Out-Patient Surgical Centers, Surgical Centers, Healthcare Facilities, ICU Areas, Isolation Wards, Emergency Rooms, Operating Rooms/Theaters and X-ray CAT Labs.

This product cleans, disinfects and deodorizes hard, non-porous hospital/medical surfaces in one step with no rinsing required.

SPORICIDAL PERFORMANCE AGAINST C. DIFFICILE SPORES: This product kills and/or inactivates the following spore in 2 minutes at 4 oz. per gal. of 400 ppm hard water, on hard, non-porous surfaces. *Clostridium difficile*.

DISINFECTANT PERFORMANCE: This product kills the following bacteria in 2 minutes at 4 oz. per gal. of 400 ppm hard water and 5% soil, on hard, non-porous surfaces: *Staphylococcus aureus* Community Acquired Methicillin Resistant (CA-MRSA) (Genotype USA300), *Escherichia coli*, *Escherichia coli* (Extended Spectrum β -Lactamase) (ESBL) (B44-156), *Klebsiella pneumoniae*, *Klebsiella pneumoniae* Carbapenem Resistant (BA4-176), *Staphylococcus aureus* Methicillin Resistant (MRSA), *Pseudomonas aeruginosa*, *Salmonella enterica*, *Staphylococcus aureus*, *Staphylococcus aureus* Vancomycin Intermediate Resistant (VISA) (HP 5836), *Enterococcus faecalis* Vancomycin Resistant (VRE), *Streptococcus pneumoniae*.

TUBERCULOCIDAL PERFORMANCE: This product kills the following mycobacteria in 10 minutes at 4 oz. per gal. of 400 ppm hard water and 5% soil at 21° C, on hard, non-porous surfaces: *Mycobacterium bovis* (Tb).

VRUCIDAL* PERFORMANCE: This product kills the following viruses in 2 minutes at 4 oz. per gal. of 400 ppm hard water and 5% soil, on hard, non-porous surfaces: Adenovirus Type 5 (Strain Adenovid 75), Hepatitis B Virus# (HBV), Hepatitis C Virus# (HCV), Herpes Simplex Type 1 Virus, Herpes Simplex Type 2 Virus, Human Immunodeficiency Virus Type 1 (HIV-1) (HTLV-IIIb), Influenza A Virus, Norovirus, Respiratory Syncytial (RSV) Virus, Rhinovirus Type 37, Rotavirus (Strain WA), Vaccinia Virus.

*Indicates a 5-minute contact time is required for this claim.

CONTINUED FROM LEFT PANEL

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product can be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

FOR USE AS A SPORICIDE AGAINST *Clostridium difficile* SPORES:

1. Pre-clean surfaces. Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the disinfectant product. This cleaning may be accomplished with any cleaning solution, including this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

2. Apply use solution of 4 oz. of this product per gal. of water to disinfect hard, non-porous surfaces with a sponge, brush, cloth, mop, by immersion, auto scrubber, mechanical spray device, hand pump trigger spray device, coarse trigger spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.

3. Treated surfaces must remain wet for 2 minutes.

4. Wipe dry with a clean cloth or allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.

5. Prepare a fresh solution daily or when visibly dirty.

NOTE: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

FOR USE AS A ONE-STEP, GENERAL, HOSPITAL, MEDICAL DISINFECTANT, VIRUCIDE*, DEODORIZER, CLEANER:

1. Pre-clean heavily soiled areas.

2. Apply use solution of 4 oz. of this product per gal. of water to disinfect hard, non-porous surfaces with a sponge, brush, cloth, mop, by immersion, auto scrubber, mechanical spray device, hand pump trigger spray device, coarse trigger spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.

3. Treated surfaces must remain wet for 2 minutes.

4. Wipe dry with a clean cloth or allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.

5. Prepare a fresh solution daily or when visibly dirty.

FOR USE AS A TUBERCULOCIDE:

1. Pre-clean heavily soiled areas.

2. Apply use solution of 4 oz. of this product per gal. of water at 21°C to disinfect hard, non-porous surfaces with a sponge, brush, cloth, mop, by immersion, auto scrubber, mechanical spray device, coarse trigger spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.

3. Treated surfaces must remain wet for 10 minutes.

4. Wipe dry or allow to air dry.

5. Prepare a fresh solution daily or when visibly dirty.

***KILLS HIV-1, HBV AND HCV ON PRE-CLEANED HARD, NON-POROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV AND HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.
Personal Protection: Wear protective latex gloves, gowns, masks and eye protection.

Cleaning Procedure: Blood and other body fluids containing HIV-1, HBV and HCV must be thoroughly cleaned from hard, non-porous surfaces and objects before application of this product.

Disposal of Infectious Materials: Blood and other body fluids, cleaning materials and clothing must be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

Contact Time: ≠ Allow solution to remain wet for 2 minutes for HIV-1, 5 minutes for HBV and HCV at 4 oz. of this product per gal. of water.

AUTOMATED DILUTION SYSTEM (NAVIGATOR) INSTRUCTIONS: Remove cap and insert cartridge into dispenser. See dispenser instructions for proper placement of cartridge. Press button to dispense 4 oz. per gal. of water of this product into a bucket, bottle, or other container.

Product No.
48027



CONTINUED FROM LEFT PANEL

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER, CORROSIVE. Causes irreversible eye damage and skin burns. Harmful if swallowed. May be fatal if inhaled. Do not get into eyes, on skin or on clothing. Do not breathe vapors or spray mist. Wear goggles or face shield and chemical-resistant gloves and protective clothing when handling. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

ENVIRONMENTAL HAZARDS

This pesticide is toxic to birds, fish and aquatic invertebrates.

PHYSICAL OR CHEMICAL HAZARDS

STRONG OXIDIZING AGENT, CORROSIVE. Mix only with potable water at 60 - 80° F. Product must be diluted in accordance with label directions prior to use. This product is not combustible; however, at temperatures exceeding 150° F, decomposition occurs releasing oxygen. The oxygen release could initiate combustion. Never bring this product into contact with other sanitizers, cleaners or organic substances.

FIRST AID

In case of emergency, call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

IF ON SKIN: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15 - 20 minutes.

IF SWALLOWED: Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person. Call a poison control center or doctor immediately for treatment advice.

IF INHALED: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth, if possible.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

PESTICIDE STORAGE: Store only in original container. Keep this product under lock and storage sufficient to make it inaccessible to children or persons unfamiliar with its proper use.

PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance.

CONTAINER HANDLING: Non-Refillable Container. Do not reuse or refill this container. Wrap empty container and put in trash or offer for recycling.

(J-17)

Manufactured for:

SS TRIPLE S
2 EXECUTIVE PARK DRIVE,
BILLERICA, MA 01862

CONTINUED ON RIGHT PANEL