



PERISEPT

Sporicidal Disinfectant Cleaner

Two-Minute
C. diff Claim
at 99.9999%
Efficacy



New Power Against COVID-19, C. diff and other Common Pathogens

SSS Perisept Sporicidal Disinfectant Cleaner kills SARS-CoV-2, the virus that causes COVID-19, in 2 minutes and also eliminates and prevents the spread of *C. difficile*, as well as, MRSA and many other common pathogens in healthcare, including ESKAPE pathogens – Enterococcus faecalis (VRE), Staphylococcus aureus, Klebsiella pneumoniae, Acinetobacter baumannii, Pseudomonas aeruginosa, and Enterobacter species.*

Contains no alcohol or bleach, requires no rinsing and leaves no film on hard surfaces. As a hospital-grade disinfectant cleaner and sporicide it is ideal for surgical rooms, patient rooms, operating suites, physical therapy departments, nursing services, autopsy facilities, long-term care facilities and nursing homes.


According to the Centers for Disease Control and Prevention, EPA-registered disinfectants with a sporicidal claim have been used successfully with environmental surface disinfection in patient care areas.

**When used in accordance with product labeling.*



EPA Registration
10324-214-12120

How does PERISEPT Work?



Combat *C. diff* and HAIs with the Industry's Fastest Sporidical Disinfectant Cleaner

Fast Kill Time

With its shorter kill times, SSS Perisept Sporidical Disinfectant Cleaner stays wet and works faster. The outcome – assured efficacy against hard-to-kill *C. difficile* spores and other dangerous pathogens every time you clean and disinfect.

Cost-Effective Non-Bleach Formula; Superior Material and Surface Compatibility

The non-bleach, hydrogen peroxide/peroxyacetic acid based formulation leaves no harmful or corrosive residues, no undesirable film on surfaces. SSS Perisept Sporidical Disinfectant is compatible with many different surfaces and medical equipment used in healthcare including stainless steel, aluminum, chrome, glazed ceramic tile, plastic and painted surfaces, and finished floors. Unlike bleach-based products, it will not corrode metal or damage mattress covers, and requires no rinsing. See the Technical Guide for a broad list of compatible surfaces.

Benefits of a Concentrate

SSS Perisept Sporidical Disinfectant is designed for use with the closed loop Navigator Dilution Control System. The Navigator DCS offers closed loop dispensing and minimizes contact with the concentrate. The easy to follow numbering and color-coded labeling simplifies training for EVS departments. The high 1:32 dilution rate offers customers a lower end-use cost making it affordable for use as a daily cleaner, in addition to operating rooms, isolation rooms and post discharge. Standardizing to one disinfectant cleaner reduces training challenges, potential confusion, as well as, the number of chemicals stored in the supply closet.

PERISEPT Dispensing Options



#23200

Navigator 4 - Station Dispenser

Four-Station Wall Mount Dilution Control System accurately dispenses four different concentrated cleaning products in one easy step.

Features:

- Easy refill bottle installation
- Single dial product selector
- One-handed spray bottle filling
- Remote bucket/auto scrubber fill gun with trigger on handle
- Rinse water selector
- 2.5 GPM dispensing with action-gap back-flow protection
- No metering tip installation—refill bottles come pre-tipped
- Locking doors prevents access to concentrates

Safe-Guard™ System All Navigator Models

Leak-proof bottle insert integrates with dispenser connector—eliminates spills; prevents contact with concentrated product.



#23215

MPD Wall Mount

Compact design wall-mounted bucket fill and bottle fill dispenser

- Multiple Product Dispenser (MPD) - easily switch between a broad range of Navigator Concentrate Refills
- One hand fill for spray bottles and remote fill gun for bucket fill
- Fill spray bottles, auto-scrubbers, mop buckets, microfiber buckets and other containers



Bucket Fill
#23202

Bottle Fill
#23201

Single Station Units

Efficient, accurate dilutions in an ultra-compact design

- Available in bottle fill or bucket fill models
- Durable polypropylene construction with optional locking cabinet
- Safe-Guard System eliminates spills and chemical contact
- 2.5 GPM dispensing with action-gap back-flow protection
- Dispensers can be connected to a single water source



#23203

Navigator Portable

Saves time—no back-and-forth trips to refill buckets or equipment

- Ideal for remote locations – works with any standard water source
- Compact, rugged design with simple one-hand control
- Safe-Guard System eliminates spills and chemical contact

PERISEPT

Sporicidal Disinfectant Cleaner



ITEM#	DESCRIPTION	SIZE	YIELD/CASE	PACK
48027	Perisept Sporocidal Disinfectant Cleaner	2 Liters	68 quarts	2/cs
SSS Navigator Dilution Control System Solutions				
23200	Wall Mount-Four Station Dispenser, includes water supply hose			1/cs
23201	Wall Mount-Single Station Dispenser, Bottle fill*			1/cs
23202	Wall Mount-Single Station Dispenser, Bucket fill*			1/cs
23215	MPD Wall Mount Compact Dispenser, Bottle fill & Bucket fill*			1/cs
23203	Portable Dispenser*			1/cs
23209	Premium Grade Water Supply Hose, 6 ft., Black			1/cs
SSS Peritex EVS Wipe Solutions				
81004	Peritex EVS Wipes, 10"x12", 1/4 Fold, White, 320 ct.		40/pk, 8 pks/cs	
81005	Peritex EVS Wiping System Bucket			6/cs
81006	Peritex EVS Wipes, 10"x12", Rolls, White, 320 ct.		40/roll, 8 rolls/cs	
81007	Peritex EVS Wiping System Canister			6/cs
SSS NexGen Microfiber Solutions				
19090	Blue Microfiber Cloth, 16" x 16"			24/cs
19091	Pink Microfiber Cloth, 16" x 16"			24/cs
19092	Green Microfiber Cloth, 16" x 16"			24/cs
19093	Yellow Microfiber Cloth, 16" x 16"			24/cs
19100	PreTreat Orange Bucket for Clean Cloths, 3 gal.			3/cs
19101	PreTreat Gray Bucket for Soiled Cloths, 3 gal.			3/cs
19108	Microfiber Cleaning Dolly II			1/cs
33001	Healthcare Cleaning Cart for EVS (JCAHO Compliant)			1/cs

*Requires Water Supply Hose



Distributed by: